**Truking's Overall  
Solution for Freeze-Dried Preparation**  
楚天冻干制剂整体解决方案

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Equip global pharmaceutical manufacturing industry with intelligent wings.

纯化水制备装置 Purified Watergenerator Equipment	- 05 -	<b>无菌药品冻干制剂生产整体解决方案</b> - 01 - Overall Solution for Sterile ProductsLyophilization Production
注射用水制备装置 WFI Generationdevice	- 06 -	
纯蒸汽发生器 Pure SteamGenerator	- 07 -	
水分配系统 Water DistributionSystem	- 08 -	
		<b>制药用水系统</b> - 04 - Pharmaceutical Water System
		<b>全自动配液系统</b> - 09 - Fully Automatic Preparation System
脉动真空蒸汽灭菌柜 Fully Automatic Preparation System	- 11 -	<b>灭菌系统</b> - 11 - Sterilizing System
器具清洗机 Appliance Cleaning Machine	- 14 -	<b>胶塞（铝盖）清洗机</b> - 18 - Rubber Stopper(Cap) Washing Machine
订制化清洗架 Customized Cleaning Rack	- 16 -	
立式胶塞（铝盖）清洗机 Vertical Type Rubber Stopper(Cap) Washing Machine	- 18 -	
卧式胶塞（铝盖）清洗机 Horizontal Type Rubber Stopper(Cap) Washing Machine	- 20 -	<b>西林瓶洗烘灌封联动线</b> - 21 - Vial Compact Line for Washing, Sterilizing & Drying, Filling & Stoppering, Capping
立式超声波清洗机 Vertical Ultrasonic Washing Machine	- 22 -	<b>自动进出料系统</b> - 25 - Auto Loading & Unloading System
隧道式灭菌干燥机 Sterilizing & Drying Tunnel	- 23 -	
灌装加塞机 Filling & Stoppering Machine	- 24 -	<b>真空冷冻干燥机</b> - 33 - Vacuum Freeze-Dryer
固定式进出料系统 Fixed Feeding and Discharging System	- 25 -	<b>冻干实验室</b> - 35 - Lyophilization Laboratory
阿尔法（α-ALUS）系列自动进出料系统 (α-ALUS)Series Automatic Feeding and Discharging System	- 27 -	<b>轧盖机</b> - 37 - Capping Machine
移动式 RGV 进出料系统 Mobile RGV Feeding and Discharging System	- 29 -	<b>全自动灯检机</b> - 39 - Fully Automatic Inspection Machine
无推杆进出料系统 No-pushrod Loading and Unloading System	- 31 -	<b>机器人自动包装生产线</b> - 41 - Robot Automatic After-Packaging Compact Line
		<b>生产制造与质量保障</b> - 42 - Manufacturing and Quality Assurance

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Founded in 2000, Truking Technology Limited has become a leading pharmaceutical equipment company in China and renowned worldwide. The company specializes in pharmaceutical equipment and integrated solutions and takes the initiative to push forward the R&D of smart pharmaceutical factories. It is a listed A-share company with multiple wholly-owned or holding subsidiaries, such as Germany-based Romaco Group, Truking Watertown, Sichuan Pharmaceuticals Design Institute, Truking Feiyun, Truking Ingenuity, Truking Micro-Sphere, Truking SUT, Truking Gene, Wachine, Truking Changxing, Truking Scientific Instrument, Truking Jingbang Engineering Technology, Truking Pitide Biotechnology, Truking Boyuan Intelligent Technology, Truking New Material, Truking Intelligent Robot etc. Truking has over 9000 employees around the globe, with a total asset of over 11 billion RMB.

Truking has two operational headquarters in Changsha, China, and Germany, with research and development institutions including Changsha Central Technology Research Institute, Europe Technology Research Institute, and Sichuan Pharmaceuticals Design Institute. It has established many platforms for technology innovation, such as the National Enterprise Technology Center, National Innovation Center, Post-Doctoral Scientific Research Station, and Academician Expert Workstation. By December 31, 2023, a total of 5,466 patent applications from China have been submitted, with 4,313 patents granted and 3,090 patents currently in force. Additionally, 54 PCT international patent applications were filed, resulting in 26 patents granted in multiple countries including the United States, Russia, India, South Korea, Germany, Indonesia, Japan, and Europe. The company has led the development of over 20 national industry product technical standards for this system. The products and services of the group have now reached over 180 countries and regions, with its international market share steadily increasing year by year.

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### Overall Solution for Sterile ProductsLyophilization Production

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
Adopt latest improved process, high preparation efficiency, low running cost and modular design. It is easy for installation, operation and maintenance.

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Introduce advanced technology from abroad, adopt new design concept; the whole process of equipment adopts full computer control. It is safe and reliable, with thoroughly sterilization effect.

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
Design as per specific requirement of customer, adopt advanced and flexible modular process piping integrated unit to realize seamless connection with filling.

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Adopt method of single arm hanging inner cylinder rotating structure, with 3 cycles of washing with SIP. There is no dead angle during washing with high washing efficiency. Auto closed unloading system, sterile transfer and conveying system can be provided as optional according to user's production mode. It has strong operability, with less manual intervention and high sterile security level

5

**洗灌灌封联动线**  
Compact Line For Washing, Sterilizing & Drying,  
Filling & Stoppering, Capping

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Compact line for washing, sterilizing & drying, filling & stoppering, capping: adopts L shape in-wall structure, one-side operation, effectively realizes function division, reduces cross contamination, improves aseptic security level. ORABS/CRABS, on-line IPC weighing system and filling CIP/SIP can be equipped as option to realize isolation and automation of auxiliary operation

6

**冻干机及自动进出料系统**  
Freeze Dryer And Auto Loading  
& Unloading System

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Freeze dryer and auto loading & unloading system: realize aseptic auto transferring from filling & stoppering machine to freeze dryer, freeze dryer to capping machine completely, reduce persons in aseptic room, reduce contamination risk, improve aseptic security level.

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Provide machine model with different function according to different cleanliness class; can realize aluminum particles auto collection, imaging testing before capping and other functions

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adopts machine to substitute person, improves inspection efficiency and reliability, has small space and low running cost, can realize automatic inspection for visible particle, appearance defect and capping defect in freeze drying products.

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effectively connect production and store, auto storing and retrieving, improve labor productivity, reduce labor

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adopts automatic packaging compact line of pallet making, pallet entering, three-dimensional wrapping, packing, strapping, palletizing This line can reduce labor intensity and production cost, improve production efficiency and packing quality rate.

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
realize production data collection of each system, person working state collection, informatization of equipment operation & maintenance management, can communicate with global customers through remote service measuring and control center, realize preventive maintenance, remote video for trouble confirmation and maintenance guidance

# 制药用水系统 Pharmaceutical Water System

## ▶ 纯化水制备装置 Purified Watergenerator Equipment

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The preparation system is designed according to the customer's individual requirements, and has an advanced concept team and a professional manufacturing, on-site construction, and verification team, meets requirement of latest version of GMP, EUGMP, FDA and cGMP. It can provide systems for following two fields

## 生物技术领域 / Biotechnology Field

- ⊙ 上游培养基配液及离心均质设备辅助工艺系统集成
  - ⊙ 下游分离纯化设备辅助工艺系统集成
  - ⊙ 血液制品生产工艺系统集成
  - ⊙ 疫苗制品生产工艺系统集成
  - ⊙ CIP/SIP
- + Upstream medium preparation and centrifugal homogeneity equipment auxiliary system integration
  - + Downstream separation and purification equipment auxiliary system integration
  - + Blood products production system integration
  - + Vaccine products production system integration
  - + CIP/SIP

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
Advanced control system, can optimize effectiveness of producing process, guarantee stability of product quality, reduce production cost, increase production flexibility. We have experienced local electrical automation department to provide customized high price performance ratio automatic control system as your special requirement, including:

- + Electrical design
- + Control cabinet manufacturing
- + FDS development
- + Software programming
- + System test
- + User train
- + System commissioning

## 灭菌系统 Sterilizing System

### 脉动真空蒸汽灭菌柜 Pulsating Vacuum Steam Sterilization Cabinet

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Pulsating vacuum steam sterilization cabinet is a machine integrating vacuum technology and fluid technology. It is composed of box, door, frame, outer sheet metal, vacuum system, steam inlet system, exhaust system, compressed air system, transfer system, filter online sterilization system, hot air system, control system and other functional branch systems. The pulsating vacuum steam sterilization process can effectively remove the non condensable gas in the chamber, so as to ensure that the steam can penetrate into the porous loose load quickly and effectively, and ensure that the temperature in the cabinet is uniform. After sterilization, the product is dried by vacuumizing and adding clean and dry air to remove water vapor in the cavity. Process flow: preheating - pulsating - heating - sterilization - drying and cooling

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

**面向厂房定制化灭菌设备 / Customized Sterilization Equipment for Factory**

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Based on the sterilization objects and quality requirements of customers, provide relevant validation services

- + Vacuum leakage test
- + The Bowie Dick test
- + No load heat distribution
- + Full load heat distribution
- + Heat penetration test
- + Drying degree test

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Appliance cleaning machine is a process equipment which integrates appliance cleaning and drying. In each cleaning step, through the coordination of cleaning temperature, pressure and time, and according to the requirements of key quality parameters, different cleaning agents are used. The cleaning process is pressurized by the clean booster pump, and 100% coverage cleaning is carried out by the spray ball in the cleaning chamber. The verified cleaning process uses low conductivity as the daily monitoring data. In the drying step, the fan is used to take the air, and the heater and high efficiency filter are used to blow and dry the inside of the pipeline, the cleaning chamber and the inside and outside surface of the appliance. Process flow: pre washing - Cleaning - rinsing - rinsing - drying and cooling.

**主要清洗对象 / Main Cleaning Objects**

不锈钢器具  
Stainless steel appliances

软管  
Hose

塑料容器  
Plastic container

玻璃器皿  
Glassware

工具  
tool

- 配置低电导率监测，对终洗水进行日常监控，确保清洗效果
  - 配置酸碱罐及自动输送系统，清洗阶段自动混入清洗剂，并配置高电导率检测，可实现清洗剂配比溶度精确控制及记录
  - 独立双路热风干燥系统，风机自动取风经加热器及各级过滤器对管路、清洗舱及器具内外表面吹扫烘干
  - 配置工控机，具备审计追踪、数据报表功能
  - 采用西门子 S7-1200 系列 PLC，运行速度快，可扩展的模块多
  - 满足 FDA、cGMP 要求，控制系统设计符合 Part 11 和 GAMP5 要求
- Equipped with low conductivity monitoring, daily monitoring of final washing water to ensure the cleaning effect
  - Equipped with acid-base tank and automatic conveying system, the cleaning agent can be automatically mixed in the cleaning step, and the high conductivity detection can realize the accurate control and record of the ratio and solubility of the cleaning agent
  - Independent two-way hot air drying system, the fan automatically takes the air, through the heater and all levels of filters to purge and dry the internal and external surfaces of pipelines, cleaning chambers and appliances
  - The industrial computer is equipped with audit tracking and data reporting functions
  - It adopts Siemens s7-1200 series PLC with fast running speed and many expandable modules
  - It meets the requirements of FDA and cGMP, and the control system design meets the requirements of Part 11 and GAMP5

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

◎ The cleaning rack is customized according to the cleaning equipment, and the cleaning coverage rate is 100%

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○ Provide relevant verification services based on customer's equipment and cleaning requirements

## 胶塞（铝盖）清洗机 Rubber Stopper(Cap) Washing Machine

### ▶▶ 立式胶塞（铝盖）清洗机 Vertical Type Rubber Stopper(Cap) Washing Machine


#### 优势特点 / Advantages and Characteristics

- 国内首创基于流体动力学原理的可控脉动悬浮清洗系统；
  - 利用真空脱泡空化效应及水汽悬浮动力对胶塞进行清洗，清洗效果优良；
  - 独立清洗罐设计，内筒与罐体合二为一，无清洗死角；罐体整体转动，内无机械搅拌，不损伤胶塞，无额外颗粒产生；
  - 罐体进出口充气密封，静态密封，无摩擦颗粒产生，泄漏风险低；
  - 腔体有效容积 85%，清洗效率高，能耗更低；
  - 多种出料方案，满足各种无菌转运需求。
  - 可处理两种物料，通过配置胶塞清洗罐和铝盖清洗罐
- + The first domestic controllable pulsating suspension cleaning system based on the principle of fluid dynamics;
  - + The rubber plug is cleaned by vacuum deaeration cavitation effect and water vapor suspension power, and the cleaning effect is excellent;
  - + Independent cleaning tank design, the inner cylinder and tank body are integrated, without cleaning dead angle; the tank body rotates as a whole, without mechanical stirring, without damage to the rubber plug, and without extra particles;
  - + The inlet and outlet of the tank are gasified and sealed statically, without friction particles and low leakage risk;
  - + The effective volume of the cavity is 85%, with high cleaning efficiency and low energy consumption;
  - + A variety of discharging schemes to meet the needs of aseptic transportation.
  - + It can handle two kinds of materials by configuring rubber plug cleaning tank and aluminum cover cleaning tank

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The mechanical arm lifts the rubber plug cleaning tank and butts with the rubber plug feeding station, as shown in the figure below:

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## 卧式胶塞（铝盖）清洗机 Horizontal Type Rubber Stopper(Cap) Washing Machine

### 优势特点 / Advantages and Characteristics

- ◎ 国内首创基于流体动力学原理的可控脉动悬浮清洗系统；
  - ◎ 利用真空脱泡空化效应及水汽悬浮动力对胶塞进行清洗，清洗效果优良；
  - ◎ 独立清洗罐设计，内筒与罐体合二为一，无清洗死角；罐体整体转动，内无机械搅拌，不损伤胶塞，无额外颗粒产生；
  - ◎ 罐体进出口充气密封，静态密封，无摩擦颗粒产生，泄漏风险低；
  - ◎ 腔体有效容积 85%，清洗效率高，能耗更低；
  - ◎ 多种出料方案，满足各种无菌转运需求。
  - ◎ 可处理两种物料，通过配置胶塞清洗罐和铝盖清洗罐
- + The triple cleaning system is designed to ensure the cleaning effect, the primary cleaning cycle cleaning, the saving of cleaning water, and the fine washing and rinsing;
  - + Quantitative silicone oil filling process makes the silicification effect more uniform;
  - + The fixed discharge valve has no risk of friction particles and secondary pollution, and can be cleaned and sterilized online;
  - + The door of the chamber is designed to be opened backward, and the maintenance work is carried out in the equipment room; the sealing ring of the door is statically sealed, and the risk of leakage is low;
  - + The vacuum pump cooling system, with fast cooling speed and adjustable capacity, can reuse the cooling water and spray cooling the hot drainage steam;
  - + Design of vacuum breaking filter module for online sterilization and filter integrity detection;
  - + Backflow prevention device for steam water separation at outlet: air partition and water seal structure are adopted, with steam pressure relief outlet to prevent pollution risk caused by floor drain backwash;
  - + The temperature and pressure are monitored in real time, and the curves are printed and recorded in real time to form a complete electronic data report;
  - + The portal automatic control system fully meets the requirements of GAMP5 and 21cfrpart11;
  - + A variety of discharge schemes to meet the needs of aseptic transportation.

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### Vial Compact Line for Washing, Sterilizing & Drying, Filling & Stoppering, Capping

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Mainly used for the production of vials or freeze-dried powder injections in pharmaceutical factories. During the linkage production, it can automatically complete more than 20 processes such as bottle cleaning, preheating, drying and sterilization, depyrogenation, cooling, filling, and stoppering.

## ▶▶ 立式超声波清洗机 Vertical Ultrasonic Washing Machine

### 优势特点 / Advantages and Characteristics

- ◎ 机架整体加工成型
- ◎ 金属烤漆处理
- ◎ 集中供油系统
- ◎ 提供管路材质证明、焊接记录及内窥镜检查报告
- ◎ 可增加注射用水降级水清洗工艺
- ◎ 选配压缩空气排空系统
- + Integral processing of the frame
- + Metal paint treatment
- + Centralized oil supply system
- + Provide piping material certification, welding records and endoscopy inspection report
- + Can increase WFI degraded process
- + Optional compressed air evacuation system

## ▶ 隧道式灭菌干燥机 Sterilizing & Drying Tunnel

### 优势特点 / Advantages and Characteristics

- ◎ 多路 PID 温控技术
  - ◎ 分风框及模块化回风框结构
  - ◎ 负压密封
  - ◎ 风压自动平衡系统
  - ◎ 停电系统
  - ◎ 网带 CIP 欠载
  - ◎ 灭菌模式
  - ◎ 冷却段冷却方式：风冷 or 水冷
- + Multi-channel PID temperature control technology
  - + Air distribution frame and modular air recycling structure
  - + Negative pressure sealing
  - + Automatic wind pressure balance system
  - + Power-off protection
  - + CIP for the belt
  - + Sterilization mode
  - + Cooling method of cooling section: air cooling or water cooling

## ▶▶ 灌装加塞机 Filling & Stoppering Machine

### 优势特点 / Advantages and Characteristics

- ◎ 采用直线或 L 型设计
  - ◎ 全伺服控制，确保运行稳定性和精度
  - ◎ 无瓶不灌、无瓶不加塞
  - ◎ 自动剔废
  - ◎ 选配在线称重系统
  - ◎ 选配 CIP/SIP
  - ◎ CRABS/ISOLATOR
- + Adopt straight line or L-shaped design
  - + Full servo control to ensure running stability and accuracy
  - + No bottle no filling, no bottle no stoppering
  - + Automatic rejection
  - + Optional online weighing system
  - + Optional CIP / SIP
  - + CRABS / ISOLATOR

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The feeding and discharging of freeze-drying machine is a kind of fully automatic linking equipment which uses automation technology to complete the semi stoppered vials from filling machine to freeze-drying machine, and from freeze-drying machine to capping machine on the premise of meeting the regulations. It has the advantages of high mobility, good repeatability, easy verification, stable and reliable product quality, and minimizing product pollution. According to different implementation forms, it can be divided into fixed automatic feeding and discharging system and mobile automatic feeding and discharging system.

## 优势特点 / Advantages and Characteristics

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
双排进料 (550pcs)  
Double row feeding (550pcs)

双排蛇形缓冲  
Double row snake buffer

波纹管在线完整线检测  
On line complete line inspection of bellows

星形轮过载保护  
Star wheel overload protection

Crabs 腔体穿越伸缩网带  
Crabs cavity passing through telescopic mesh belt

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Alpha ALUs series automatic feeding and discharging system is a solution for the future. It not only has the characteristics of small space occupation, flexible layout, good aseptic compliance, intelligent and efficient feeding, but also can be combined with lyotk industrial, pilot, small-scale freeze-drying machines and different Rabs forms to form an overall freeze-drying solution. It can effectively reduce the input cost of the equipment and maximize the use value of the whole line.

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### Flexible layout solution

- The area of grade B area is reduced by 70%, and the operation cost is reduced
- Open flow of people and logistics channels
- It can meet the feeding demand of single or multiple lyophilizers
- It is suitable for lyotk industrial freeze-drying machine, small-scale and pilot scale freeze-drying overall solution
- It can integrate orabs, crabs and Isolator systems

### Aseptic compliance

- Embedded feed push rod module
- Design of bellows without push rod
- Clean control of key areas

### Efficient and stable feeding mode

- Feeding speed: 600vmp
- Continuous pushing mode
- Double push design, single row pushing and whole plate feeding work division and cooperation
- Big star wheel design to improve the efficiency of bottle counting

## ▶▶ 移动式RGV进出料系统 Mobile RGV Feeding and Discharging System

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Mobile RGV feeding and discharging system adopts RGV technology to separate bottle collecting station, transfer station and discharging station. Through RGV car to complete the freeze-drying machine feeding and discharging action. The system includes the whole line feeding system, mobile loading and unloading trolley system, mobile track and power supply system, discharging system, Rabs laminar flow system, etc;

## 优势特点 / Advantages and Characteristics

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RGV 小车  
RGV car

出料平台  
Discharge Platform

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## No-pushrod Loading and Unloading System

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Mainly used for automatic loading and unloading of lyophilized APIs in trays, especially suitable for high-end biological lyophilized APIs and other products with small batch, high added value, and biological activity. The RGV technology is combined with the tray lifting structure to realize automatic loading without push rods.

## 优势特点 / Advantages and Characteristics

- 结构简单，便于清洁；
  - 占地面积小，操作简单；
  - 无前推波纹管，无菌风险小；
  - 配合托盘输送机构，可实现全自动化进出料。
- + Simple structure, easy to clean;
  - + Small size, easy cooperation;
  - + Design of bellows without push rod, sterile and low risk;
  - + Tray conveying mechanism, realizes fully automatic loading and unloading.

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Mainly used in the fields of pharmaceuticals, food and health care products such as biological products, chemicals, bottle freeze drying (lyophilized injection), plate freeze drying (raw material), natural medicine, heat-sensitive medicine, antibiotics, oral freeze-dried tablets, realizes freeze drying process.

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
**Chamber:** Formed at one time to reduce the risk of leakage. The chamber adopts square structure, and all processed surfaces and holes are formed by advanced processing equipment at one time, reaching the current international advanced level; Adopts large rounded corner structure of R> 50mm; the single-angle tilt of the bottom is more favorable for drainage; 316L material for all inner material; inner surface polished to mirror surface <0.4umRa.

**Product shelves:** The shelf is made of AISI 316L stainless steel; Adopts vacuum brazing to ensure good flatness and sealing of the shelf; 316L stainless steel hose with braided mesh is used; the surface is polished to 0.4umRa and passivated, and adopts the helium mass spectrometer leak detection to ensure no leakage

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**Condenser:** The Condenser can be designed as a vertical type or a horizontal type, also it can be placed on the side/back/bottom of the chamber; Mushroom valve design for best vapour flow and ice distribute; 316L material for all inner material, inner surface polished to <0.6umRa. Condenser coil surface area design based on 1.1:1 with the shelf area.

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**Chamber door:** Special hinge design to ensure the Main door move easily, and equipped with D-form gasket for good sealing

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Process development: Help user to develop freeze drying curve

Recipe design: Help user to develop recipe of freeze drying product

Process Scale up: Offer help to user who need to transform from small scale

production or laboratory to large scale production

Process optimization: Help user who has original process curve, but need to optimize, to optimize the process curve, improve freeze drying speed, and reduce freeze drying cycle.

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
Mainly used for capping of vials in pharmaceutical factory. It adopts side-placing singleroller sealing structure, can complete the following process automatically: vial infeed, cap infeed, capping. As per requirement of new version GMP, capping can be made in A grade ventilation condition under B, C or D grade background. When capping in C or D grade background, it needs auto detection of sealing for stopper and vial mouth before capping.

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Mainly used for automatic inspection of visible particle, apparent defect and sealing defect for vial freeze drying products in pharmaceutical factory.


## 产品特点 / Product Features

- ◎ 产品缺陷: 冻干块上表面、下表面以及侧表面的异物, 如玻璃、纤维、黑块、装量过高、过低或者空瓶冻干块萎缩、液化、粉化, 冻干块歪斜或倾倒。
  - ◎ 封口瑕疵: 有无胶塞、轧盖(铝盖)位置, 有无刮痕, 弯曲或者凹痕、有无塑料盖、塑料盖是否破损、有无黑点。
  - ◎ 瓶身瑕疵: 玻璃有无划痕、裂纹或结石瓶身上是否有污点。
- + Product defect: Foreign bodies on the upper, lower and side surfaces of the lyophilized block, such as glass, fiber, black block, filling volume too high, too low, or empty bottle, lyophilized block shrinking, liquefied, powdered, inclining or dumping.
  - + Sealing defect: Whether there is stopper, cap position, whether there are scratches, bends or dents, whether there are plastic covers, whether the plastic cover is damaged, and whether there are black spots
  - + Vial body defect: Whether there are scratches, cracks or stone on the glass, whether there is dirty spot

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Station 1: lower surface inspection of freeze-dried block

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
Station 2: side surface inspection of freeze-dried block

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Station 3: upper surface inspection of freeze-dried block

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Station 4: inspection of vial body defects, including scratch, crack etc

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Station 6(optional): static detection of caps(detection for cap top)

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Mainly used for package of ampoule injection, vial injection, powder injection, lyophilized powder injection, oral solution, large infusion, tablets, capsules and other medicines, as well as cosmetics, food and other items after filling and sealing. Can automatically complete the process of light inspection, labeling, pallet making, loading, boxing, wrapping, case packing, bundling, palletizer supervision code management and others.

# 生产制造与质量保障 Manufacturing and Quality Assurance

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

As a key public project of the pharmaceutical factory, the water system follows the V model from planning to operation, the document meets verification requirements, and the automatic control system can meet GAMP

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### + “Housekeeper” service

Rich experience on project management;  
Strong team for project management;  
Professional & Abundant talent reserve

### + Intimate pre-sale service

Overall solution design as per user' s requirement  
Consulting for overall solution of user' s factory  
Free technical training from experts

### + Professional sales service

Re-optimization design adjustment for customer required product  
Operator training, guidance  
Installation, commissioning and site training

### + Responsible after-sales service

24-hours consulting service  
Convenient and quick on-site service  
Equipment and technology upgrading service

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## Service Concept

400-9988-900

Adhering to concept of “all for customers, for all customers” , and with the mission of “creating value for customers” , TRUKING focuses on building first-class R&D, design, sales, manufacturing, service teams and service management system to offer high-efficient and high-quality products and services for the customers.

### ◎ Project Planning Support

TRUKING has a series of professional medicine and pharmacy design institutes, product research institutes, a group of senior scientific research and technical personnel, aseptic pharmaceutical process verification experts, professors, etc., to provide scientific and effective technical planning support for customers' projects.

### ◎ Process Assurance Service

TRUKING launched the overall solution for medical & pharmaceutical equipment realizing the full automation and aseptic production process, intelligent detection and sorting of the drug, as well as computer verification assurance so as to produce assured products for the market.

### ◎ Service In The Park

As an equipment supplier, TRUKING not only integrates industrial regulations and standards into product design and manufacturing, but also introduces the industrial regulations and standards into Truking industrial park in the form of lectures and exchanges to interact and integrate with the technology and production of pharmaceutical enterprises.

### ◎ Equipment Management Consultation Service

TRUKING aims to satisfy the maximum benefited demands of the customer and to provide a package of equipment management services to realize minimization investment and maximization production output together with pharmaceutical factories.

◎ Adhere to the philosophy of "customer-centric", TRUKING creates the extraordinary service experience with the spirit of the craftsman, to create more value for customers.

◎ We are always ready to escort the normal operation of the customer's equipment.

- + Onsite equipment installation, commissioning, acceptance
- + Training for equipment process, operation and maintenance
- + Periodic door-to-door maintenance and return visit
- + Performance test and maintenance for machinery parts
- + Timely supply of spare parts
- + Product technology upgrading and secondary development
- + Equipment overhaul

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